



IRMRA

INDIAN RUBBER MANUFACTURERS RESEARCH ASSOCIATION

Affiliated to Ministry of Commerce & Industry (Government of India)

(A NABL (ISO - 17025-2005) Accredited Laboratory)

E-mail : info@irmra.org • Website : www.irmra.org

EVALUATION REPORT

Report No : RPT / 486 / 17-18	25 May 2017	Page 1 of 2
-------------------------------	-------------	-------------

--- CUSTOMER DETAILS ---

Customer Name & Address :	
M/s.Vijayalakshmi Health and Surgicals P.Ltd Plot No:406, Apiic Limited, Growth Center, Gullapali, Ongole, Prakasham Dist, Andra Pradesh-523211	
Contact Person :	R Bhuvaneswari
Telephone No :	-
Customer Letter Ref No with Date :	Dated 22 May 2017
Date of Receipt :	22 May 2017
Any Other Information like witness / third party	

--- SAMPLE DETAILS ---

Sample Code / Description	Surgical Hand Gloves
Product Specification	-

--- REPORT DETAILS ---

Description of Report	Page Nos
RPT / 486 / 17-18 - M & R -	--
RPT / 486 / 17-18 - C & I -	2 of 2

Terms and Conditions

1. Sample was not drawn by IRMRA, unless otherwise stated.
2. The result enclosed are only related to the sample submitted to us.
3. The report shall not be reproduced except in full without prior written permission of laboratory.
4. Information contained in this letter / or enclosed reports is/are purely meant for your guidance and should not be used for litigation/evidence/advertisement, without prior permission from laboratory.

भारतीय रबड़ निर्माता अनुसंधान संघ
Indian Rubber Manufacturers Research Association

Prepared and Checked By
एन.ए.फोंडके (विशेष अधिकारी-ग्राहक सेवा)
N. A. Phondke (Senior Officer-Customer Service)



भारतीय रबड़ निर्माता अनुसंधान संघ
Indian Rubber Manufacturers Research Association

Suchismita Sahoo
Authorized Signatory
सुचिस्मिता साहू (सहायक निदेशक)
Suchismita Sahoo (Assistant Director)



IRMRA

INDIAN RUBBER MANUFACTURERS RESEARCH ASSOCIATION

Affiliated to Ministry of Commerce & Industry (Government of India)

(A NABL (ISO - 17025-2005) Accredited Laboratory)

E-mail : info@irmra.org • Website : www.irmra.org

TECHNICAL REPORT / C & I

RPT / 486 / 17-18 - C & I

25 MAY 2017

Page 2 of 2

Sample : Surgical Hand Gloves Specn. -

Test (Condition)	Method	Unit	Specified Value	Observed Value
Extraction (n-Hexane/Reflux temp./first 7hrs)	US-FDA-177-2600. Clause-F	mg/inch ²	175Max	0.95
Extraction (n-Hexane/Reflux temp./ Succeeding 2hrs)	US-FDA-177-2600. Clause-F	mg/inch ²	4Max	0.20

Note: Sample Passes the test.

Prepared and Checked By



भारतीय रबड़ निर्माता अनुसंधान संघ
Indian Rubber Manufacturers Research Association

मनोहर नवाले (वरिष्ठ वैज्ञानिक अधिकारी)
Manohar Nawale (Senior Scientific Officer)